Piranha/Cobia & Quick Check

Getting Started

Piranha/Cobia & Quick Check Getting Started - English - Version 1.14A

Welcome to Quick Check

Ocean 2014 with Quick Check is a is a powerful tool for everybody working with Quality Assurance of X-ray systems. Ocean 2014 can be used with the X-ray meters Piranha and Cobia from RTI.



INDEPENDENT X-RAY QUALITY ASSURANCE

NOTICE

RTI Group AB reserves all rights to make changes in the Ocean 2014 software, and the information in this document without prior notice. RTI Group AB assumes no responsibility for any errors or consequential damages that may result from the use or misinterpretation of any information contained in this document.

Copyright © 2010-2019 by RTI Group AB. All rights reserved. Content of this document may not be reproduced for any other purpose than supporting the use of the product without prior permission from RTI Group AB.

Microsoft, Microsoft Excel, Microsoft Access, Windows, Win32, Windows 95, 98, ME, NT, 2000, XP, 2003, Vista, Windows 7 and Windows 8 are either registered trademarks or trademarks of Microsoft Corporation in the United States and/or other countries.

OpenOffice.org and OpenOffice.org Calc are registred trademarks of OpenOffice.org.

BLUETOOTH is a trademark owned by Bluetooth SIG, Inc., USA.

Contact Information -World-Wide

RTI Group AB Flöjelbergsgatan 8 C SE-431 37 MÖLNDAL Sweden

Phone: Int. +46 31 7463600

E-mail Sales: sales@rtigroup.com Support: support@rtigroup.com Service: service@rtigroup.com

Web site: https://www.rtigroup.com

Contact Information -United States

RTI Electronics Inc. 33 Jacksonville Road, Building 1 Towaco, NJ 07082 USA

Phone: 800-222-7537 (Toll free) Int. +1-973-439-0242 Fax: Int. +1-973-439-0248

E-mail

Sales: sales.us@rtigroup.com Support: suppor.ust@rtigroup.com Service: service.us@rtigroup.com

Web site: https://www.rtigroup.com

Intended Use of the Ocean Software

Together with instruments from RTI Group AB the Ocean Software is intended to be used for independent service and quality control, including measurements of kerma, kerma rate, kVp, tube current, exposure time, luminance, illuminance, and dose area product, within limitations stated below.

If installed according to accompanying documents, the product is intended to be used together with all diagnostic X-ray equipment except for:

- therapeutical X-ray sources.
- X-ray equipment with tube potential below 18 kV or above 160 kV.
- X-ray equipment on w hich the instrument cannot be mounted properly.

- specific types of X-ray equipment listed in the instructions for use or in additional information from the manufacturer.

With the X-ray installation without patient present, the product is intended to be used:

- for assessing the performance of the X-ray equipment.
- for evaluation of examination techniques and procedures.
- for service and maintenance of the X-ray equipment.
- for quality control of the X-ray equipment.
- for educational purposes, authority supervision etc.

The product is intended to be used by hospital physicists, X-ray engineers, manufacturer's service teams, and other professionals with similar tasks and competencies. The operator needs training to be able to use the product as intended. This training can be achieved either by study of the manual, study of the built-in help function in measurement softw are or, on request, by a course ordered from the manufacturer.

The product is intended to be used inside X-ray rooms ready for clinical use and can safely be left sw itched on and in any measuring mode in the vicinity of patients.

The product is NOT intended to be used:

- for direct control of diagnostic X-ray equipment performance during irradiation of a patient.
- so that patients or other unqualified persons can change settings of operating parameters
- during, immediately before, or after measurements.

- for any guidance to diagnosis of patients.

Table of Contents

1	Quick Check	. 6
	1.1Start the Tablet	6
	1.2First measurement with Quick Check	8
	1.3Run Quick Check on a Tablet	12
	1.4Change Meter Settings	19
	1.4.1 Piranha	. 19
	1.4.2 Cobia	. 19
	1.5Applications with analysis	28
	1.6Quick Check Favorites	31
	1.6.1 Detector selection	. 31
	1.6.2 Create Quick Check Favorites	. 31
	1.7Save a measurement	33
	1.8Open an old measurement	34
	1.9Print	34
	1.10. Transfer data to Excel	36
	1.10.1 Excel control window	. 36
	1.10.2 Connect (Fixed format)	. 36
	1.10.3 Connect (Standard mode)	. 36
	1.10.4 Send data to Excel	. 36
	1.11.Import/export	45
	1.12. Options	47
	.1.13. Close Quick Check	50



Use Quick Check when you just need to make a quick measurement. Quick Check automatically adapts to the meter you use and the detector you have connected. It is fully plug-and-play and any change you do is automatically detected and Quick Chick adapts to the new situation. Quick Check is designed for a tablet (touch screen) but works of course also on a standard laptop. You can do the following with Quick Check.

- Quick measurements
- Save
- Print
- Export to Excel
- Applications with analysis

1.1 Start the Tablet

This is a quick introduction to how you get started with the tablet that comes with your meter. The tablet has Windows 8 and is configured and ready to be used with the meter.

To start measuring:

- 1. Power on the meter.
- 2. Power on the Tablet.



3. Windows 8 starts and the Start screen is launched.

Start			rti 🎴
Cozen 2014	Ка	SyDire Pictos	Stype
RTI Updater 2014 RTI Manuals	Celendar	Resting List	Food & Drink
	Prople	Store	News
	rope	5004	News

4. You will see the Ocean 2014 tile to the left, click on it to start Ocean 2014.



- 5. Assuming that your meter is powered on; Ocean 2014 starts in Quick Check mode.
- 6. Select the modality you want to use:

0							Ocean 201	4					-	8 ×	
Reset	Start	0 Pause	Copture	Meter settings	Position check	New Quick Check	Application	Save	Preview/Print	Export		Options	Close Quick Check	() Help	
🗸 Conn	ected														
C Heasure	ments 🃫 No i	leasurelle	nt loaded			S	elect modal	lity		×					
					➔ Radio	graphy									
					➔ Fluoro	scopy									
					Mamn	nography									
					→ CT										
					Dental										
· Piranh	a						Battery 78 9	%							
	< P										[11111]		o	4:29 PM	8

The coming topics will describe how you perform your first measurement and how you navigate in Quick Check.

How to turn off the Tablet

It is recommended that you turn off the Tablet, and not just puts it in sleep mode, when you don't use it. It consumes power even in sleep mode and the battery will be drained quite fast even if you don't use it.

To turn off the tablet:

1. Press the power button and keep it pressed until the screen shows the following:



2. Release the button and slide the panel down to shut down the Tablet.



1.2 First measurement with Quick Check

By default, Quick Check is launched automatically when Ocean 2014 starts and the meter is available. Use Quick Check when you just want "a meter" to do a quick measurement. Quick Check is extremely simple to use since it automatically recognizes what type of meter you have and which detectors you use. It is full plugand-play, directly when you unplug or plug in a detector; Quick Check sets up the measurement for you. You can save data you measure with Quick Check but you must save manually - Quick Check will never ask you if you want to save. It is also possible to print and export measured data to Excel. Quick check also has built-in applications for accuracy, linearity, reproducibility, half-value layer and more. Quick check is designed for use with a tablet and easy to operate directly on a touch screen.

Let us start with a simple measurement to illustrate how to use Quick Check. This example assumes a Piranha 657. If you have another model follow this example and apply it on what you see on your screen when you use your Piranha.

Make sure your meter is powered on and within distance for Bluetooth communication or connected via USB. Disconnect any external probe that might be connected (unless you have a Piranha model that only works with external detector or a Cobia Sense). Assume that we are using a Piranha 657.

1. Start Ocean 2014 or start the Quick Check from the Measure tab on the ribbon bar if you are in Ocean 2014 main view: (Note: If you use Windows 8; you start Ocean 2014 from the desktop)



2. Quick Check starts and a menu will be shown. The list is depending on your meter type, it might not appear at all.

elect	modality
•	Radiography
÷	Fluoroscopy
÷	Mammography
e	СТ
¢	Dental

3. Select for example Radiography. For other modalities from here and on; it may be several choices to make before the measurement screen appear. For example, for mammography you must chose calibration and if compression plate is used or not.

4. The measurement loads and you are ready to start to measure.

0	Untitled - Ocean 2014 Professional		_ 🗆 🛛
Reset Start Pause Capture	Meter settings Position deck New Quick Check Application Favorites	Save Preview / Export	Options Close Quick Help Check
1¢ Ready Measurements 1¢ Untitled			
Tube voltage 0 kV W/3 mm Al	Exposure time 0 ms W/3 mm Al	200 Tube voltage	Exposure rate
Exposure 0 mGy W/3 mm Al	Exposure rate 0 mGy/s ^{W/3 mm Al}	100 1.0 - 140 (9)(5)(1,4 - (¥))(20- 1,2 - 49 100 1.2 - 49 100	
HVL 0 mm Al	Total filtr. 0 mm Al W/3 mm Al		
		0.2 20 0 12 14 15 15 0 10 12 14 15 Time (m)	20 22 24
		Show/hide Cursor 1 Cursor 2 Diff Tube voltage Exposure rate Time 0,000 0,000 0,000	kv mGy/s ms
	I		
Piranha 45 - 125 kV	Battery 100 %	Messages	

If the selected Quick Check is using a calibration that isn't available in your meter/detector, a dialogue is shown. You can here select another calibration instead of the one saved in the document.



In some situations a message is shown with a yellow background; all such messages can be closed by just clicking or tapping on the yellow background.



5. Make some exposures.

0	Untitled - Ocean 2014 Professiona	al	- 🗆 🗙
Reset Start Pause Capture	Meter Position New Quick Check Application Favorites	Save Preview / Export O	ptions Close Quick Help
♥ Ready ↑ Untitled			
Tube voltage 87,27 kV ^{W/3 mm Al}	Exposure time 99,87 ms W/3 mm Al	25- 90- 80- 80- 80- 80- 80- 80- 80- 80- 80- 8	Exposure rate
Exposure 2,414 mGy ^{W/3 mm Al}	Exposure rate 23,68 mGy/s W/3 mm Al Total filtr.		
3,56 mm Al	3,4 mm Al		
			100 120
		Showy muce Cursor 1 Cursor 2 Diff Tube voltage 87,37 37,37 0,00 kt Exposure rate 23,59 23,48 0,1054 m Time 24,19 96,76 72,57 m	(g) Gy/s S
 Piranha 45 - 125 kV 	Battery 100 %	Messages	

6. You can operate Quick Check with mouse and keyboard or with your fingers if you have a tablet with touch screen.

Quick Check is designed to be simple to use on a tablet (touch screen) but you can of course also use a mouse. Below are the most important operations described.

Toggle between displays and logged data: Click on the button in the lower left corner of the display panel or sweep left or right on the displays if you use a tablet.

Enlarge waveform: Click on the button in the lower left corner of the waveform graph or double-click on the waveform. Same to minimize it.

Zoom waveform: Zoom the waveform by holding down the left mouse button and mark the part you want to zoom. When you have zoomed it is remembered for current exposure. If a new row is automatically created for the next exposure, the zoom state is inherited for the new row.

Enlarge a display/single display: Double-click on a display.

Change unit of measure: Right-click on a unit and choose a new unit from the list shown. To right-click one tablet; hold your finger on the unit, remove it when a small square appear.

Change kV range: Click on the kV range indication in the lower part of the screen.

You have now seen the basic function of the Quick Check. You can now try to connect an external detector (if you have one) and see that Quick Check will directly recognize the probe and ask you what to do.

Read also the full description of Quick Check in the topic Run Quick Check on a Tablet to further see how you operate Quick Check. This section is followed by more about how to store, print and how to export your measured data to Excel.

Messages

It is indicated on the lower status bar if a message arrives.

	- rose ronage
	Exposure rate
	Time
Battery 100 %	Messages

When a message arrives the icon changes to:

Messages (1 new)

You can now click on the icon to open the message list:

0	Messages	- 🗆 🗙
Date	Subject	
2015-01-02	The Calibration of your instrument S/N CB2-13060173 is Due on 2015-02-01	
Delete		Close

Here are all your messages shown. Double-click in the list to read a message. Unwanted messages can be deleted.

You can turn off calibration reminders in the Program Options, see more in section Program Options

1.3 Run Quick Check on a Tablet

Quick Check is simple to operate and use both on a tablet with touch screen or a normal laptop. This section will describe all the actions you can use to operate Quick Check if you use a tablet, as well as a laptop.

Quick Check main screen for standard measurements looks like this:

0	Untitled - Ocean 2014 Professiona	al	_ 🗆 🛛
Reset Start Pause Copture	Meter settings Position New Quick Check Application Favorites	Save Preview / Export	Options Close Quick Help
珍 Ready で Measurements 却 Untitled			
Tube voltage O kV W/3 mm Al Exposure	Exposure time O ms W/3 mm Al Exposure rate	200 Tube voltage	Exposure rate
0 mGy W/3 mm Al	0 mGy/s W/3 mm Al Total filtr.	$g_{1,0}^{1,0+}$ 40 - $g_{1,1}^{1,0+}$ $g_{120}^{1,0}$ - $g_{1,2}^{1,2+}$ $g_{100}^{1,0}$ - $g_{1,0+}^{1,0+}$ g_{100-} -	
0 mm Al	0 mm Al W/3 mm Al	808+2 80- 0.8+ 00- 0.4- 40- 0.2- 20-	
		0 1 0 12 14 16 15 Show/hide Cursor 1 Cursor 2 Diff ↑ Tube voltage ► Exposure rate	20 22 24
■ ✓ Piranha 45 - 125 kV	Battery 100 %	Tme 0,000 0,000 0,000	ms

You have a ribbon bar at the top with different buttons to activate different function. You can click on these buttons with a mouse click or by tapping on the screen with your finger if you use a tablet.

Switch between display and list with exposures

Click on the button in the lower left corner or sweep on the display area of the screen.

)	Untitled - Ocean 2014 Profession	al	
Reset Start Pause Capture	Meter Position New Quick Check Application Favorites	Save Preview / Export	Options Close Quick Help
19 Ready			
Tube voltage 87,27 kV	Exposure time 99,87 ms		Exposure rate
W/3 mm Al Exposure 2,414 mGy	Exposure rate 23,68 mGy/s	(20 - 70 - (20 -	
W/3 mm Al HVL 3,56 mm Al	W/3 mm Al Total filtr. 3,4 mm Al	816-880 50-40 50-40 610-42 40-40 610-42 30-	
W/3 mm Al	W/3 mm Al		80 100 120
Sweep left	or right to toggle	Time (mi) Time (mi) Show/hide Cursor 1 Cursor 2 Diff • Tube voltage 87,37 87,37 0,00 • Fronsure rate 23,59 23,48 0,1054	kv kv
between di	splays and log	Time 24,19 96,76 72,57	ms
Piranha 45 - 125 kV	I Battery 100 %	🐼 Messages	

13

Click on the button to toggle between displays and log

The log is shown and you can see all the measurements you have done during this Quick Check session.



Click on the button to toggle between displays and log

Sweep or click on the button to go back to the displays.

Enlarge the waveform

You can maximize the waveform by double-clicking (tap twice) on it.





Zoom the waveform

Tap on the waveform and move finger (or right-click, hold down and move mouse pointer) to mark the zoom.



15

The zoomed waveform is shown when you release and you can now see details.



The current zoom state is now automatically locked and the waveform will be saved in this way. If a new row is automatically created, the new row inherit the zoom state.

You can use the scroll bar to "move" backward and forward. Click on the button in the right upper corner to go back to full view.

Waveform data and cursors

The waveform data is shown under the display (when it is minimized). You can switch between "cursor values" and "value between cursor", use double-click (or tap twice) to change. Values for the two vertical cursors are shown. Cursor data for the horizontal cursor is only shown in the graph.



Single display

If you just want to see one value, you can maximize just that display. Double-click (tap twice) on the display you want to see.





Double-click on the display to go back to smaller displayes.

Change unit of measure

You can change unit of measure by right-clicking on the unit:





Select a new unit.

Change calibration

It is possible to change calibration by right-clicking in the display. This is not applicable for radiography, fluoroscopy and dental since for these modalities only one calibration (W/AI) is required.

D	Ur	ntitled - Ocea	n 2014 Prot	fessional					-	
Reset Start Pause Capture	Meter settings	New Quick Check	Application	Favorites	s Save	Preview / Print	Excel	Options	Close Quick Check	? Help
🍄 Ready										
Tube voltage O kV Selenia/MCP-Mo/Mo Exposure O mR Selenia/MCP-Mo/Mo HVL O mm Al Selenia/MCP-Mo/Mo Bight-click on the	Exposure time 0 Selenia/MCP-Mo/Mc Exposure rate 0 Selenia/MCP-Mo/Mc	ms R/min		Ennosume rate (Rinin)		Tube voltage		3	bposure rate	
To right-click on the tablet: Hold down your finger on the unit until a square is vis mathen remove it.				le,	ið ihow/hide Tube voltage Exposure rate Time	12 14 Curso	16 13 Time (m r1 Cursor	s 20 s) 2 Diff 0.000	22 2 kV R/min s	4 (1) ~

For mammography selection of calibration may look i two different ways; either as a list of manufacturer and models with specific calibrations or as a list of all available calibrations. The list with manufacturer may look like this when you start a new Quick Check:

Select manufacturer		
➔ Fujifilm	Select model ×	Select beam quality
GE Healthcare	🔿 Selenia	Selenia/MCP-Mo/Mo
→ Hologic	Selenia Dimensions	➔ Selenia-Mo/Rh
IMS Giotto	MultiCare Platinum	→ Selenia-W/Rh
Siemens	Affirm Prone	→ Selenia-W/Ag
Generic	All calibrations	➔ Show all
All calibrations	🔿 Back	🔿 Back
→ Back		

The list with all calibrations may look like this (the exact look depends on actual calibrations available in the meter):



Change number of decimals

You can change number of decimals in the same way by right-clicking on the numeric value in a display.

1.4 Change Meter Settings

Piranha and Cobia have slightly different meter settings but the way you access and change them is the same. Cobia has less settings, since it has auto-range and is a less advanced meter than Piranha. Quick Check will recognize which meter you use and what its capabilities are and adapt to it. You will only see what you can use and you can without any changes switch between Piranha and Cobia.

W/3 mm Al	W/3 mm Al
4,76 mm Al	3,0 mm Al
W/3 mm Al	W/3 mm Al
Piranha 99 - 155 kV	💷 Charging
	5 5

Currently meter is shown in the lower left corner. Quick Check can be directly used with both Piranha and Cobia, just connect to the meter you want to use.

Note: You can use the same Ocean 2014 templates for both Piranha and Cobia, a template built for Piranha can be used with Cobia and vice verse. It is recommended, if you intend to use your templates with both Piranha and Cobia, that you build your templates for Piranha. The reason for this is that Piranha has more settings and you can in this way setup the templates to work in the best way with Piranha. f you do the opposite (you build it for Cobia), Ocean 2014 will when you use it with Piranha select default settings for the Piranha when a set value is missing (due to it doesn't exist for the Cobia). If you build templates without a meter connected and intend to use it with both Piranha and Cobia; go to program options and select default "Meter type" in the Preference section.

1.4.1 Piranha

There are several ways to change meter settings. All meter settings are available if you click on the "Meter adjust" button:

Ocean 2014			_		Conc. Concerna	Acres 1	100 mil 1-4
Reset	Start	Pause	Capture	Meter setting	s Position check	New Quick Check	Application
🏟 Ready							
💼 Measureme	nts 🧳 NEV	/ Quick-Check 4	exposures				
Tube welte						-	
Tube volta	ge				Exposure tim	ie	
1	20,02	2 kV			1	00,4 n	ns
W/3 mm Al					W/3 mm Al		

When you click on the button the Meter settings are shown (this screen may look different depending on model and selected detectors):

0	Untitled - Ocean 2014 Professional – 🗆 🗙					
Reset Start Pause	Capture Meter settings	Image: Constraint of the deck Im				
Pready						
Internal detector	kV sensitivity	General Delay (ms) Window (ms) Post delay (ms)				
45 - 125 kV ~	High ~					
Calibration W/3 mm Al	Dose/TF sensitivity Low v	Measuring mode Waveform recording time Normal v 1024 ms (2 kSa/s) v				
Beam correction factor		Waveform type				
External detector						
Not used	Not used Trigger level (time) 50 %					
□ Normalized exposure and exposure rate SDD: cm SSD: cm						
→ right sweep to close						
✔ Piranha 45 - 125 kV		Charging al				

There are also some short cuts to change some of the set values. You can click (or tap) on the kV range indication at the bottom to change kV range:

W/3 mm Al	W/3 mm Al
	lotal filtr.
4,76 mm Al	3,0 mm Al
W/3 mm Al	W/3 mm Al
Piranha 90 - 155 kV	💷 Charging

When you click or tap here, a menu is shown that allows you to quickly select a different kV range.

When you measure on mammography you might want to change calibration and/or status for the compression paddle:

0		
Reset Start Pause Capture Meter	settings Position check New Quick Application Save Preview,	Print Export
12 Ready		
Measurements of Ontilled		
Tube voltage	Exposure time	
Tube voltage		110-
0 kV	0 ms	108-
Mo/20 um Mo	Mo/20 um Mo	106-
Exposure	Exposure rate	
		104-
0 mGy	0 mGy/s	102-
Mo/30 um Mo	Mo/30 µm Mo	100-
HVL		
		98-
0 mm Al		96 -
Mo/30 µm Mo		94 -
	— \	
		92-
Click or tap here to	Right-click calibration	90 -
	text in any display to	9,0 9,2 9,4
change status for	text in any display to	
compression naddle	change calibration	Show/hide Cur
compression paddie	× v	I ube voltage
		Exposure rate
		Time
Piranha 18 - 49 kV Compression Paddle	e used 🖾 Charging	

The following meter settings are available in Quick Check:

Meter setting	Description and use
Delay	Add a delay after the detection of trig before measurement of kVp starts. This will delay the kVp measurement, it doesn't affect dose, mAs or time measurements.
Window	If a time is specified, kVp is measured during the window time (starts after the delay)
Post delay	This is the time the meter waits after trig off before it assumes that the exposure is finished. The post delay must be set to a time longer than any dead time in the radiation.
kV range	Current kV range. You need to change this for radiography/fluoroscopy and CT, mammography and dental have only one range.
Calibration	Available calibrations for the internal (kVp and exposure) used. See table below that describes usage of the different calibrations.
kV sensitivity	Sensitivity setting for the kVp detector (internal detector). Hi = High sensitivity - for low dose rate Lo = Low sensitivity - for high dose rate
Dose/TF sensitivity	Sensitivity setting for the dose and total filtration measurement (internal detector). Hi = High sensitivity - for low dose rate Lo = Low sensitivity - for high dose rate
Sensitivity (External)	Sensitivity setting for the external probe Hi = High sensitivity - for low signals Lo = Low sensitivity - for high signals
Beam Correction factor	General (user-defined) correction factor used for all exposure related parameters measured with the internal detector.
Beam Correction factor (External)	General (user-defined) correction factor used for all exposure related parameters measured with the external detector.

Meter setting	Description and use
Calibration (External)	Calibration for the external probe
Measuring mode	You can select between Normal, Timed or Free run Normal =use this measuring mode for exposures and fluoroscopy Timed = meter measures during a time you specify Free run = meter measures continuously without use of any trig levels (for more information see table below)
Measurement time (Timed mode only)	Measuring time when Timed mode is used.
Waveform recording time	Select the waveform recording time. Use the shortest time to see details in the waveform. If you use a longer time, you lose details in the waveform. This setting doesn't influence on the accuracy.
Waveform type	This is the waveform type for the X-ray generator. It is normally HF/DC. Sometimes for older X-ray units and for dental you must use 1-phase. Be careful to select the correct waveform type for maximum accuracy. Note that there is a special selection for AMX-4.
Total filtration (External) (R/F, dental only)	The filtration used the external RTI Dose Probe to do energy compensation.
Added filtration (mammography only)	Added filtration used for the internal detector do do energy compensation and kV compensation.
Added filtration (External) (mammography only)	Added filtration used for the external RTI Dose Probe do do energy compensation.
Compression paddle (mammography only)	Select here if compression paddle is used or not used.
Compression paddle thickness(mm) (mammography only)	Specifies the thickness of the compression paddle in equivalent mm Al. Default is 0.12 mm.
Trigger level (time)	This is the level used for the time measurement. You can use this if you want for example to avoid pre-pulses to be included in the exposure time.
Normalized exposure and exposure rate	The square-law is used to normalize the dose and dose rate to a distance (SSD) different from the measuring distance (SDD) when the box is checked.

How to use the different calibrations:

Code	Calibration	Usage	
R1 C1	W/3 mm Al	General radiography, fluoroscopy, dental and CT	
C3	Straton (Siem1)	Suitable for Siemens CT with Straton tube	
C4	GECT (7°)	Suitable for GE CT tubes with a 7° anode angle as well as for other manufactures CT tubes and replacement tubes with a 7° anode angle	
C5	Aquillion 64-	Suitable for Toshiba Aquillion 64-320 CT	
C6	GECT (10.5°)	Suitable for GE CT tubes with a 10.5° anode angle	
C7	GECT (Cardiographe)	Suitable for GE CT Cardiographe	
M1	Mo/30 µm Mo	General mammography	
M3	Mo/25 µm Rh	General mammography	
M4	Rh/25 µm Rh	General mammography	
M6	W/50 µm Rh	General mammography - suitable for Hologic Selenia Dimensions and Fujifilm Amulet	
M7	W/0.50 mm Al	General mammography - suitable for Philips MicroDose (Sectra)	
M8	Mo/1 mm Al	General mammography	
M10	M10 W/50 μm Ag General mammography - suitable for Hologic Selenia Dimensions Fujifilm Amulet		
M11	W/75 µm Ag	General mammography	
M12	W/50 µm Rh (Gio)	Suitable for Giotto Mammography	
M15	5 W/0.70 mm Al General mammography - suitable for Hologic Selenia Dimensions a Fujifilm Amulet		
M16	W/50 µm Ag (Sel)	Suitable for Hologic Selenia	
M17	W/50 µm Rh (Sel)	Suitable for Hologic Selenia	
M18	W/0.30 mm Cu	General mammography - Suitable for Hologic Selenia Dimensions and Fujifilm Innovality/Cristalle	
M19	W/0.70 mm Al (Inno/ Crist)	Suitable for Fujifilm Innovality/Cristalle	
M20	W/50 µm Rh (Inno/ Crist)	Suitable for Fujifilm Innovality/Cristalle	
M21	Mo/25 µm Rh (Sel)	Suitable for Hologic Selenia	
M22	Rh/30 µm Ag (GE HC)	Suitable for GE Senographe Prestina	
M23	Rh/30 µm Ag IQST (GE HC)	Suitable for GE Senographe Prestina	
M24	Mo/0.25 mm Cu (GE HC)	Suitable for GE Senographe Prestina	
M25	Rh/0.25 mm Cu (GE HC)	Suitable for GE Senographe Prestina	
M26	Mo/30 µm Mo (GE HC)	Suitable for GE Senographe Prestina	
M27	Affirm Prone W/Ag	Suitable for Hologic Affirm Prone	
M28	Affirm Prone W/Al	Suitable for Hologic Affirm Prone	

There are three different measuring modes available using the Piranha. They are as follows:

Measuring mode	Description and use
Normal	The Normal mode is used for short and long (fluoro) exposures. In this mode, your meter will automatically sense if there is a signal and when it is above a certain trigger level. If the exposure is long, the displays/grid will be updated with new data every 2 seconds. If the exposure is short, the results are displayed as soon as the trigger is off.
Free run	The free run mode has no trigger level. As soon as the meter is told to begin measuring, it starts to measure even if there is no signal. This measuring mode is useful when the signal you want to measure is very low. Free run is recommended for light measurements, especially when measuring "ambient" light (when no shutter is present).
Timed	The Timed mode setting measures during a pre-defined time period. Measurements in Timed mode must be started manually. This measuring mode is very useful when you want to measure a very low signal. You can use the "very high" sensitivity setting in Timed mode and it will further improve the meter's capability to measure very low signals.

You can read more about meter settings for Piranha in the topics Meter Adjust and Different measuring modes.

1.4.2 Cobia

There are several ways to change meter settings. All meter settings are available if you click on the "Meter adjust" button:

Ocean 2014				0	Concernance of the local division of the loc	Annual State	and and in such
Reset	D Start	Pause	Capture	Meter settings	Position check	New Quick Check	Application
🗘 Ready	IP Ready Image: Measurements Image: Measurements Image: Measurements Image: Measurements						
Tube volta	ge			I	Exposure tim	e	
120,02 kV				100,4 ms			
W/3 mm Al					W/3 mm Al		

When you click on the button the Meter settings are shown:

•		Untitled - Ocean 2014 Pro	fessional			_ 0 ×
Reset Start Pause C	pture Meter Position check	New Quick Check Application	Favorites	Preview / Excel	Options	Close Quick Help
Pready						
Internal detector (Cobia CB4-140 Calibration Beam cor W/3 mm Al	rection factor No	ternal detector t used				
General Measuring mode Delay (ms) Win. (Normal Auto Auto Waveform recording time 7 s (3000 Sa/s) Normalized exposure and exposure	ms) Post delay (ms) 250 ms - Al Trigger level (time) 50% - e rate	MX4 waveform type				
SDD: inches SSD:	inches					
a Cobia 40-150 kV		@ 0 %	🐼 Messa	qes		

There are also some short cuts to change some of the set values. You can click (or tap) on the kV range indication at the bottom to change kV range:

	W/3 mm Al	W/3 mm Al	
	HVL	Total filtr.	
	4,76 mm Al	3,0	0 mm Al
	W/3 mm Al	W/3 mm Al	
L			
	Piranha 90 - 155 kV		Charging

When you click or tap here, a menu is shown that allows you to quickly select a different kV range.

When you measure on mammography you might want to change calibration and/or status for the compression paddle:



The following meter settings are available in Quick Check:

Meter setting	Description and use
Delay	Add a delay after the detection of trig before measurement of kVp starts. This will delay the kVp measurement, it doesn't affect dose, mAs or time measurements.
Window	If a time is specified, kVp is measured during the window time (starts after the delay)
Post delay	This is the time the meter waits after trig off before it assumes that the exposure is finished. The post delay must be set to a time longer than any dead time in the radiation.
Calibration	Calibration for the internal probe.
Beam Correction factor	General (user-defined) correction factor used for all exposure related parameters measured with the internal detector.
Beam Correction factor (External)	General (user-defined) correction factor used for all exposure related parameters measured with the external detector.
Calibration (External)	Calibration for the external probe.
Measuring mode	You can select between Normal or Timed Normal =use this measuring mode for exposures and fluoroscopy Timed = meter measures during a time you specify (for more information see table below)
Measurement time (Timed mode only)	Measuring time when Timed mode is used.
AMX-4 Waveform type	Check this box when you measure on a GE AMX-4.
Total filtration (External) (R/F, dental only)	The filtration used the external RTI Dose Probe to do energy compensation.
Trigger level (time)	This is the level used for the time measurement. You can use this if you want for example to avoid pre-pulses to be included in the exposure time.
Normalized exposure and exposure rate	The square-law is used to normalize the dose and dose rate to a distance (SSD) different from the measuring distance (SDD) when the box is checked.

How to use the different calibrations (only one calibration is available for Cobia):

Code	Calibration	Usage
R1	W/3 mm Al	General radiography, fluoroscopy and dental

There are two different measuring modes available using the Cobia. They are as follows:

Measuring mode	Description and use
Normal	The Normal mode is used for short and long (fluoro) exposures. In this mode, your meter will automatically sense if there is a signal and when it is above a certain trigger level. If the exposure is long, the displays/grid will be updated with new data every 2 seconds. If the exposure is short, the results are displayed as soon as the trigger is off.
Timed	The Timed mode setting measures during a pre-defined time period. Measurements in Timed mode must be started manually. This measuring mode is very useful when you want to measure a very low signal. You can use the "very high" sensitivity setting in Timed mode and it will further improve the meter's capability to measure very low signals.

You can read more about meter settings for Cobia in the topics Meter Adjust and Different measuring modes.

1.5 Applications with analysis

There are also a number of applications available in Quick Check. An application is a predefined measurement where your data are analyzed, it might also include pass/fail criteria. The applications available in Quick Check are not possible to modify in any way. If you have more needs, use Ocean Professional. You can then create tests with user-defined calculations and pass/fail criteria. You will find applications for:

- Accuracy
- Reproducibility
- mA linearity
- HVL
- CTDI

1. You start an application by clicking on the Application button (if the button is not enabled, no applications are available for current mode):



2. A menu showing available applications is shown:



3. Chose an application, for example kVp Accuracy. The screen layout is different for the application compared to a standard Quick Check measurement. In this case there are no displays only grid (exposure list) is shown. There is also an analysis that perform some kind of calculation (in this case kVp accuracy) on your measured data.



4. You can modify set values in the grid (exposure list). Click in a cell to change set values. For example, you may want to change the kVp set value for exposure #4 from 80 to 81 kV.

🚹 Measureme	ents	🦻 Untitled		
View / Select	#	Set kV (kV)	Tube voltage (kV)	kVp diff %
<u>ь</u> 🕨	1	50		
Ŀ.	2	60		
In.	3	70		
Ŀ.	4	80		
In.	5	90		
Ŀ.	6	100		
In	7	110		
۰				

5. Now make the exposures.

0					Untitled - Ocea	n 2014 Profes	sional	al	- 🗆 🗙
Reset	Start	Pause	Capture Se	Meter Position check	New Quick Check	Application Fav	orites	Save Preview / Export Options	Dose Quick Help Check
Ready	nts 📫 Untitled								
View / Select h, h, h, h, h, h, h,	# Set kV (kV) 1 50 2 60 3 70 4 80 5 90 6 100 7 110	Tube voltage (kV) 47,98 59,33 71,31 80,87	kVp diff % -4,0 -1,1 1,9 1,1	Exposure (mGy)	Exposure time (ms) 99,88 99,86 niz98,86	to		Econus rate Tube vo 80 + 22 70 + 10 80 + 22 80 + 22	
Provide a state of the state of							00 120		
Diff from set value (4	• • •	• • • •		• <u>1</u> &p				Show/hide Cursor 1 Cursor 2 Diff ● Exposure rate 20,14 20,000 mGy/s ● Tube voltage 81,47 80,86 0,61 kV Time 24,39 97,56 73,17 ms	(1)
🥏 Piranha	45 - 125 kV				Battery 100)%	ß	🐼 Messages	

You can maximize/minimize the exposure list section and the analysis section with double-click or by clicking on the buttons ().

If there is both displays, grid and analysis present the button is shown. Click here to toggle between displays and the grid.



The analysis part will now show the result.



6. Click on Preview/Print to print your report.

This analysis is fixed and you can not change its structure. If you have more needs and want to customize your analysis, change pass/fail criteria, use your own calculations and so on; then use Ocean 2014 Professional. The topic Analysis (Definitions) describes which analysis that are available in Ocean 2014.

You can open any Quick Check measurement as a Real-time display and in that way modify and extend its capability. See topic Open an old measurement.

1.6 Quick Check Favorites

The Favorites function in Quick Check is an optional way for you to customize the use of Quick Check. It allows you to use and quickly access templates that performs specific measurements. You use Ocean 2014's template design functions to create your own templates and you store them in a specific folder to make them available in Quick Check.

By default, no Favorites are available in Quick Check. You can create your own Favorites to be used in Quick Check by using Ocean 2014's template design capability. The section Create Quick Check Favorites describes how this is done.

If you have stored one or more favorites in the Favorite folder (in the Library) the Favorite list is shown when you start Quick Check. The list shows the Favorites you have stored in the folder.

Se	lect fa	vorite	23
	→ 1	Input Doserate	
	ب	kVp Accuracy (Mammography)	
	ب	More	

In this case are two favorites stored in the Favorites folder. If you select one, it loads and you can start to measure directly. If you click on "More...", the normal Quick Check setup procedure starts.

You can go to the Favorite list whenever you are using the Quick Check by clicking on the Favorites button.



The Favorite list will be shown and you can quickly select a another favorite.

1.6.1 Detector selection

When you load a Favorite the required detectors are checked when it loads. If a detector specified, or a compatible one, isn't available in the template, the Detector Selection window is shown. It will shown available detectors. It will be directly shown if you disconnect/connect a detector.

etector selection	n mana l	23
Parameter	Detector	
Tube voltage	Piranha 👻	
Exposure	Piranha 🔹	
P	Piranha	
	Citer -	
Status		

Click to minimize

Select the appropriate detector and click OK.

You have two choices if you for some reason isn't able to select a detector:

Do not use: In this case is the corresponding measured parameter(s) skipped.

Keyboard: In this case is a dialogue shown after each exposure that allows you to enter a value manually.

1.6.2 Create Quick Check Favorites

Perform the following steps to use Favorites in Quick Check:

1. Quick Check uses Real-time Display templates. You need to use Ocean 2014's template design function to create templates for you favorite measurements. This is described in the Design section.

2. Your templates can be ordinary "displays" or have analysis/calculations (you must have Professional to create this type) like the Quick Check Applications. The following design rules apply for Quick Check:

- To create a "display" favorite: design the Real-time display template with displays and no analysis.

- To create a favorite that is shown as an "Application": design the Real-time display template without displays and include an analysis.

- The template title you specify is shown in the Favorites list.

3. Copy or store the favorite templates to the The Favorites folder. It is located in the Library.



You can, if you have many Favorites of different types, create sub folders. Name the sub folders and store your favorites according to the way you have chosen to categorize them. The sub folders name will be shown in the Favorite list.

1.7 Save a measurement

Quick Check will not ask you if you want to save your data, you must yourself decide if you want to save. Click on the Save button in case you want to save your measured data.

Ocean 2014								0	- franciska i	
.=	D	0	0	•	- 111	0	Ð			
Reset	Start	Pause	Capture	Meter settings	Position check	New Quick Check	Application	Save	Preview/Print	
🏟 Ready										
💼 Measurem	ients 🍄 NEW	/ Quick-Check 4	exposures							
Tube volt	age			E	xposure tim	e				1
1	120,02	2 kV			10	00,4 n	าร		1	8

When you click on the Save button you are asked to enter a name for your data. Your data will be saved in the "Quick Check" folder in the Measurement database.

🖺 Measurements 🎝 Untitled	
E	
Search Q C*	
 ✓ "■ Quick Checks □ □ Inbox □ □ 2013-09-20 15:04:27 □ □ □ Check of kVp 	

Once you have saved your data the first time, Quick Check will ask you before it discards any data.

1.8 Open an old measurement

You can re-open a Quick Check measurement that you have saved. You can open it if you want to continue to measure or if you want to print it. From Ocean 2014 you can also open a Quick Check measurement in offline mode (no meter connected) to review or print the measurement.

To open a Quick Check measurement:

1. Click on Measurement tab.

Measurements 🏟 Check Of kVp	
E	
Search Q 💎	
▲ -	Measurement tab

2. Go to the "Quick Checks" folder, double-click on the measurement you want to re-open.

You can also open a Quick Check measurement as an Real-time display in Ocean 2014 if you want to modify it. To open as a real-time display:

- 1. Open the measurement folder.
- 2. Right-click on the Quick Check measurement you want to open as a real-time display.

Ocean 2014	1000			
.=		0	0	•
Reset	Start	Pause	Capture	Meter setting:
🤹 Ready				
🚹 Measuren	ients 🎝 Unt	itled		
🚞 Quick Cheo	ks			
⊿ = Quick Ch	necks 3-09-20 15:04:27	,		
		Open as Real-tir	ne display	\geq
	. [Delete		
	S	elect all		
	F	Preview		

3. The measurement will open in Ocean 2014 as a real-time display. You can now modify the measurement, what you can do depends on the license level you have. If you open a Quick Check measurement that has analysis with only Connect license, the analysis will not be visible.

1.9 Print

You can print your measured data (or create a PDF file). By default no waveform are included in the report. To include waveforms do the following:

1. Go to the list view.

n Ocean 2014								
. <u>*</u>	0				→ ¹ / ₂	<u>1</u>	0	
Reset	St	art Pau	se Captur	e Meter	settings Positio	n check	New Quick Check	Арр
🏟 Ready								
C Measurem	ents	🧳 Check Of kV	'P					
View / Select	#	Tube voltage (kV)	Exposure time (ms)	Exposure (mGy)	Exposure rate (mGy/s)	HVL (mm Al)	Total filtr. (mm Al)	
(🗠)	1	59,56	99,86	0,4184	4,086	2,28	3,0	
	2	79,32	99,89	0,7575	7,360	3,03	3,0	
ER .	3	89,95	99,86	0,9524	9,256	3,43	2,9	
1 ²⁸	4	120,02	100,4	1,552	15,01	4,76	3,0	
L.	5							
			Cli	ck o	n the y	wav	eform	h
								÷.,
			N WI	nen 1	the "R	" IS \	visble	э;
			Ìwa	vefo	rm is	inclu	uded	
			in	the r	oport			
				ule i	eport.			

2. Now click on the Preview/Print button.

New Quick Check	Application	Save	Preview/Print	Export	

3. A preview is now shown.

Ocean 2014 Report										23
3 ♣ ♣ ② 90% • ○ □ I4 4 1 ► ►	Close									
	Print date: 2013-1	1-15						Central Designment		ŕ
	Radiography									
	Test date: 2013-11-	4								
	Measurements									Ξ
	# Tubevoltage (kV)	Exposure time (ms)	Exposure (mGy)	Exposure rate (mGy/s)	HVL (mm Al)	Total filt: (mm Al)				
	1 59,58	99,95	0,4184	4,088	2.28	3,0				
	3 89,95	99,96	0,9524	9,258	3,43	2,9				
	Record #1	2,4	Equiper to	<u>ter</u>	5/0 F2	tu: olap.	Equire	•		ļ
		and and a second a second a second a second a	******				~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~~			
		d d d d The(m)	å å .å .	о 	8/2 E4		al do al al - Tradest	the state of the s		
Page 1 of 1	1:1		**************************************			has align				

4. You can now just preview the report and return, save it as a PDF file or print it.



Report template

You may want to change the report layout, for example, use you own logo. You can make you own report layout and assign it as default for Quick Check. You can do this in Program Options and you can find instructions in section Create a report template.

1.10 Transfer data to Excel

You can transfer data from Quick Check to Excel in three different ways:

- Connect to an Excel workbook. Measured data, are for each exposure, transferred to Excel starting in active cell. Used template determine data format.
- Connect to an Excel workbook. Measured data, are for each exposure, transferred to Excel starting in active cell. Use a fixed data format (compatible with Excel templates used with Xi and X2 View).
- Send current measured data to Excel. The data dump starts in active Excel cell.

Connect (Standard mode)

When this mode is used, data are transferred to Excel in the same order as the columns are appearing in the template. Only numerical values are sent to Excel, no units or waveforms.

	(¥ + ↓			-	
Ar	kiv Sta	rt Infog	a Sidlay	out For	mler Dat	a Gr
Klis	Nip Klip Kop tra ✓ ✓ Häm	p ut iera ▼ nta format	Calibri F <i>K</i> U	• 11 •	т А́ А́ <u>⊘</u> т <u>А</u> т	= ;
	Urklipp	Es.		Tecken	L.	i l
	B8	•	. (=	f_{x}		
	А	В	С	D	E	F
1						
2						
3						
4						
5		2,380663	12,28536	189,7393		
6		2,004946	12,28946	160,0896		
7		2,397109	12,33938	191,7014		
-						

Data are send according to how columns appear in the template.

Connect (Fixed format mode)

This is same as "Standard mode" but data is sent in a different format. The format is compatible with the format used with Xi and X2 View from RaySafe. Same Excel templates can be used with none or minimal modifications. This mode is only available in Quick Check. It is activated by checking "Fixed format (Excel connect)" on the Option page (in Quick Check). The format is described in the section Connect (Fixed format).

E	ج ا							
Arl	kiv Sta	art Info	ga Sidla	yout Fo	ormler D	lata Gra	inska \	∕isa
Klist	Klip E Kop ra V Här	p ut iera ▼ nta format	Calibri FKU	+ 1 ! + ⊞ +	1 - A /	• ≡ ≡	■ ≫	, → ≣
	Urklipp	Fa		Tecken		F2		J
A1		•	× v	f_{x}				
	А	В	С	D	Е	F	G	
1								
2								
3								
4								
5		8,532853	mGy	66,92717	mGy/s	0,114927	S	
6		11,99338	mGy	75,38516	mGy/s	0,146527	s	
7		12,57796	mGy	75,69736	mGy/s	0,154096	s	_

<u>Send</u>

In this mode the entire content in the template is dumped to Excel.

X	Bok1 -	Microsoft Excel		
Arkiv Start Infoga Sidlayout Formler Data Granska Visa			Ocean 2014	× 3
Calibri <u>11</u> A A = = = >>-	Radbryt text	Allmänt -	0 - 0 0 10 0 + 0	
Klistra Kopiera · FKU · Si · A · E E E E E	-a- Centrera över kolumner -	- % , * 0 c Villkorssty		
in • • Hamta format	ring 5	Tal S	🍄 Ready	
R29 • (* fx				*
A B C D E F G	H I J	K L M	N O P Q R	S T U
1 Measurements				<u> </u>
2				
3 Tube voltz Exposure Exposure Exposure HVL(mm / Total filtr.				
5 79.3241 99.88505 0.757539 7.359927 3.03317 3.026695				
6 89,95454 99,85621 0,952392 9,255628 3,425962 2,937374				
7 120,0175 100,3631 1,551991 15,00877 4,758177 3,002295				
8				
10 Waveforms				
11 Exposure #1 Exposure #2	Exposure #3			
12 Tube softage Exposure rate Tube softage Exposure rate	Tube votage	Exposure rate		
13 grand and a start and a sta	\$" 5 ⁸⁰			
	E - 20-			
	240-			
17 8, - 5, - 5, - 5, - 5, - 5, - 5, - 5, -	8 2 - ⁴ 20-			
18 of of the second sec				
19 Time (ms) Time (ms)		Time (ms)		
21 Exposure #4 Exposure #5				
22 Tube voltage Exposure rate 200 Tube voltage Exposure rate				
23 8 15 5100				
24 g & m g 15 + 8 m				
		1 4		
Klar				100% - +

Read the coming sections for detailed information about the different Excel modes.

1.10.1 Excel control window

When data is transferred to Excel a small Excel control window is shown to simplify the interaction between Quick Check and Excel. The exact content in the window is depending on mode (send or connect) and/or used detectors.

	Cap	oture	
	Start	Pause	
F Toggle between Quick Check and Excel	Cick on Start b External detector W/3 + 20 m Sensitivity Excel sheet Blad1 A2 ↓	Veasuring time is 10 seconds) utton to start measuring. m Al Low Get active cell	Measured values are shown here. These data are useful during fluoroscopy exposures since no intermediate data are shown in Excel, only the result after the exposure terminates. Look here to see when data are stable.

The most common meter settings are available (depending on used detector(s).

To close the Excel connection, click "x" in the upper right corner of the Excel control window. Both the window and the used Excel workbook will be closed. If necessary, you are asked to save.

1.10.2 Connect (Fixed format)

Data is sent to Excel using a fixed format that is compatible with the format used in Excel workbooks used with Xi and X2 view. The format is described in detail below.

To connect to the workbook:

First go to Options in Quick Check and make sure that the checkbox "Fixed format (Excel connect)" is checked.

- 1. Activate a Quick Check display.
- 2. Click on the Excel button.

Save Preview/Print Excel

3. A dialogue is shown, select "Connect to workbook".



4. Next you have to chose a workbook. You can chose a new workbook or one you already have on your computer. In this case select "New workbook".



Excel starts and an empty workbook is loaded.

5 ·	<i>∂</i> - +																										
kiv St	tart Info	ga Sidl	ayout P	ormler [Data Gri	anska Vi	sa Utv	vecklare	Tillägg	PDF Archited	t 4 Creator	♀ Berā	tta vad du vi	ll göra													
Klip	pp ut	Calibri	×	11 × A°.	A* = =	= ** -	📑 Ra	idbryt text		Allmänt	v			N 1	ormal	Bra		Dålig	Ne	utral	Anteck	ning	Ocean 201	4			
tra 💞 Hãi	mta format	F K	u • 🖽 •	- 🖄 - 🗚	• = =		🗏 🛱 Ce	entrera över	kolumner 🔹	S • %	• % 🎎	Villkors	styrd Form	abell -	räkning	Förklar	ande	Indata	Ko	ntrollcell	Länkad	cell 🚽	0.4		0 #	×	Sök oc market
Urklipp	p na		Tecken		6		Justering			s 13	al r						Forma	at					0 kV 0 m	Gy 0 mGy/s	0 ms		
	• : :	x v	f _x																				🍁 Read	(
А	В	с	D	E	F	G	н	1.1	1	к	L	м	N	0	P	Q	R	S	т	U	v	w					AC
																							Interna	detector			
																							45 - 1	25 kV		Ý	
																							W/3 r	nm Al		-	á 📖
																									High	~	<i>i</i> —
																							kV sens	livily	Law		
																							Dose/T	sensitivity	LOW	Ý	
																							Excel sh	eet			
																							Blad1			× 🦲	
																							A1	-	Get a	ctive cell	
																								_			
	-																										

5. Activate the cell, for example B5, where you want measured data to start. You can do that by entering B5 into the cell field in the Excel control box or select cell B5 in the Excel spreadsheet and click on "Get active cell".

6. Make a couple of exposures.

E									
A	rkiv	Start	Infog	ja Sidlaj	yout Fo	rmler [)ata Gra	inska V	/isa
Klis in	tra	Klipp Kopie Hämt	ut ra ▼ a format	Calibri F K U	- 11 	A .	▲ = = • = =		, → <u>=</u>
A1	Ur	rklipp •	5 :)	×	Tecken f _x		F ₂		I
	А		в	с	D	Е	F	G	
1									
2									
3									
4									
5			8,532853	mGy	66,92717	mGy/s	0,114927	s	
6			11,99338	mGy	75,38516	mGy/s	0,146527	s	
7			12,57796	mGy	75,69736	mGy/s	0,154096	S	

7. Measured data are transferred to Excel for each exposure. You can use the small minimize Excel window to control your measurement procedure while having Excel on the screen. You can change basic meter settings and the left most button is used to switch between Quick Check and Excel. Data are arriving in a fixed format (compatible with type format used by Xi and X2 View). The fixed format is defined below.

Definition of the fixed format

Multi-detector (internal detector)

There are four different detector combinations that each has its own fixed format:

- Multi-detector (internal detector) with an optional external detector
- Only an external radiation detector (Dose Probe, CTDP, T20, Ion chamber of any type)
- Only the Light detector
- Only a mAs probe

Tube voltage unit Dose unit

unit
Dose
unit
Dose rate
unit
Exp. time
unit
Pulses
unit
Empty
Empty
Pulse rate
unit
Dose/pulse
unit
HVL
unit
Tube mAs
unit
Tube mA
unit
TF
unit
Dose (ext.)

unit
Dose rate (ext.)
unit

Units are defined by Defult units in Program options in Ocean 2014 (not available in Quick Check).

Only external radiation detector

Dose
unit
Dose rate
unit
Exp. time
unit
Pulses
unit
Empty
Empty
Pulse rate
unit
Dose/pulse
unit

Units are defined by Defult units in Program options in Ocean 2014 (not available in Quick Check).

Only light detector

Lumunance or Illuminance	
unit	

Units are defined by Defult units in Program options in Ocean 2014 (not available in Quick Check).

Only mAs probe

Tube mAs			
unit			
Tube mA			
unit			
Exp. time			
unit			

Units are defined by Defult units in Program options in Ocean 2014 (not available in Quick Check).

1.10.3 Connect (Standard mode)

Data is sent to Excel using a fixed format that is compatible with the format used in Excel workbooks used with Xi and X2 view. The format is described in detail below.

To connect to the workbook:

First go to Options in Quick Check and make sure that the checkbox "Fixed format (Excel connect)" is unchecked. This ensures that standard mode is used and data transferred is defined by the actual template used in Quick Check.

1. Load a Quick Check or an Application.

2. Click on the Excel button.



3. A dialogue is shown, select "Connect to workbook".



4. Next you have to chose destination. You can chose a new workbook or one you already have on your computer. In this case select "New workbook".



Excel starts and an empty workbook is loaded.



5. Activate the cell, for example B5, where you want measured data to start. The measured values appears in the same order as the columns in the active template. No units are sent, only the numeric values. (If you want a fixed format and units; use the "Fixed format" option described in the next section).

6. Make a couple of exposures.



7. Measured data are transferred to Excel for each Exposure. The data transferred and the order of data is determined by the column present in the template used in Quick Check. You can use the minimized Excel window to control your measurement procedure while having Excel on the screen. The left most button is used to switch between Quick Check and Excel.

1.10.4 Send data to Excel

This is used when you want dump your current data to Excel. You have a choice to include or exclude waveforms (waveform are exported as pictures).

1. Assume that you have done a number of exposures and you want to dump all the measured values to Excel and that you want waveforms.

1	Ocean 2014											
	Reset	St	tart Pau	lse Captur	e Meter	settings Positio	n check	New Quick Check	Application	Save	Preview/Print	Export
	🏟 Ready				_							
	🚹 Measureme	nts	NEW Quick-	Check 4 exposure	25							
	View / Select	#	Tube voltage (kV)	Exposure time (ms)	Exposure (mGy)	Exposure rate (mGy/s)	HVL (mm Al)	Total filtr. (mm Al)			18	3+
	Ŀ.	1	59,56	99,86	0,4184	4,086	2,28	3,0				120-
	Ŀ.	2	79,32	99,89	0,7575	7,360	3,03	3,0			16	⁵ + 110-
	In.	3	89,95	99,86	0,9524	9,256	3,43	2,9			14	100-
	in. 🕨	4	120,02	100,4	1,552	15,01	4,76	3,0			6	90-1-00
	Ŀ.	5									^{™G}	2 2 80-

2. Click on the Export button. A dialogue is shown:

Select how you want to connect to the workbook
Send data to workbook Send current measured data to the workbook
 Connect to workbook Transfer data after each exposure to workbook
 Cancel Do nothing

Select "Send data to workbook"

3. A new dialogue is shown where you can chose if you want waveforms or not. In this case select "With waveforms".

×

4. Next you have to chose destination. You can chose a new workbook or one you already have on your computer. In this case select "New workbook".



5. Ocean 2014 starts Excel and opens the workbook you selected. All you measured and in this case including waveforms are transferred to Excel.



6. You will also see a minimized Ocean 2014 window that is shown. This window provides help when you work with Excel.

Ocean 2014	x
● = ≥ ● 0 ≥ ⊙ # ⊘	
17 Ready	

The left most button is used to switch between Quick Check and Excel. See topic Data link to get more information about this window. You can also read Send data to Excel to see how you can use this feature in Ocean 2014.

7. You can save your Excel workbook and further process your measured data.

1.11 Import/export

If you want to move measured data from one computer to another, an export and import function is available to you. However, only Export is available from the Quick Check mode. To use the Import function; you must go to Ocean 2014 main screen.

Export

To export one or several measurements:

1. Go to the Measurement tab.

2. Right-click on the measurement you want to export.

You can first make a multiple selection by using the Ctrl- and Shift key and then right-click to directly select more than one measurement for export.

🖺 Measurements 🦻 Untitled							
🔲 Quick Checks 🕨 🕻	Check Of kVp						
Search	ୁ ଦ୍						
Quick Checks 2013-09-20 Check Of ky	▲ ☐ Quick Checks ☐ 2013-09-20 15:04:27						
	Open Open as Real-time display						
	Delete						
	View report PDF Export						
	Select all						
	Preview						

3. Select Export. A new window is shown that allows you to here select more measurements to export. The measurements you include are shown on the right side.

Select	Select sessions and/or real	time displays	
Destination	Search P & Stee P & Real-time displays P & Quick Checks P & Inbox 2013-09-20 15:04:27 Check Of Kip	Name Check Of kVp	Type Quick Check

4. Click on Next when you have included all measurements.

Export		
Select	Destination	
Destination	Export via	
	Message	
		*
		-
		< Back Einish Cancel

6. You can chose between "File" and "E-mail". If you have an e-mail program, you can send the file directly if that is what you want. Otherwise select "File". Select a file name and destination to finish.

You can now use the file you have created to move your measurements to another computer with Ocean 2014.

Import

You can import data, but you must then close Quick Check to go to Ocean 2014 main screen and the application button.

1. Click on the Close button.



2. Ocean 2014 main screen will now appear. Click on the Application button and select Import.

(New New Save as Save as Print Export Import	Recent documents Q. Atest (Real-time display) 1. NEW Quick-Check 4 exposures (Quick Check) 2. App kVp Accuracy RAD (Quick Check) 3. kVp Accuracy (Real-time display template) 4. Check Of kVp (Quick Check)	I database F ear row Ear all Detect
	Close	Options Close Ocean 2014	

3. Locate the file you want to import and select it.

n Import	Import						
Measurements							
	Name	Туре	Destination				1
	NEW Quick-Check 4 expc	Quick Check	Inbox				
							1
		< Back	Einish		Cance	el	

4. The dialogue shows all measurements that will be imported. You will find them in the Inbox after the import is finished.



You can read more about Import and Export in the topics Import and Export.

1.12 Options

If you click on the options button a menu is shown.

\frown	L	. O X
Options	Close	? Help

You will have several options here:

Options				
Options				
Save as Favourite				
Meter info				
Make support file				
Check update				
Activate meter				
Show error reports				
About				

Select "Options".

D	Untitled - Ocean 2014 Profess	ional	- • ×
Reset Start Pause Capture Meter	Image: Second	Sove Prevew / Export Option	ns Close Quick Help Check
☆ Ready ① Heasurements ゆ Untitled			
Preferred modality • No preference (ask) Radiography & Fluoroscopy Radiography Fluoroscopy Mammography CT CT Dental	Auto. Position Check Aadiography Fluoroscopy Mammography Only for the first used target/filter For every used target/filter CT Dental	Defaults Exposure © Gray © Roentgen ✓ Auto-start with USB © Exposure assistant ✓ Show hints automatically © Confirm exit © Only use Quick Check © Enable Favorites © Always ask for Save © Always show all meters ♥ Fixed format (Excel connect)	
🖉 Piranha	🙋 100 %	a Messages	.d

There are three sections here:

Preferred modality

You can use this if you have a meter that covers many modalities but you don't want to see all choices in the Quick Check.

Auto Position Check

Specify in which situations you want Quick Check to automatically suggest a Position Check.

Defaults

Various settings that controls how Quick Check works.

Exposure unit: Select Gray or Roentgen

Language: select the language you want for the Quick Check.

Auto-start with USB: If this is checked, Ocean 2014 and Quick Check starts automatically as soon as the meter is connected with USB. Quick Check also terminates automatically when the meter is disconnected or powered off.

Exposure assistant: Turn on or off the Exposure assistant.

Show hints automatically: If a Quick Check template has a hint attached to it, the hint is shown automatically if this checkbox is checked.

Confirm exit: If you check this a dialogue is shown when you quit Ocean 2014. It give you three alternatives:

- Quit Ocean 2014 and return to Windows
- Quit Ocean 2014 and turn off computer
- Resume Ocean 2014

Only use Quick Check: If you check this box, Ocean 2014 main window is hidden and you will only see the Quick Check. Ocean only starts if meter is available.

Close the Options screen by sweeping right if you use a tablet or click on the close button.

Always ask for Save: When this box is checked Quick Check always asks if you want to save your measurements before starting a new measurement or closing.

Always show all meters: When this box is checked a list with available meters to connect to is shown. If unchecked, Ocean 2014 directly tries to connect to the last used meter.

Fixed format (Excel connect): Check this box if you want to use fixed format when connecting to Excel.

Save as Favourite...

Save current measurement as a Favourite. The name given will be shown in the list when the Favourite button is clicked.

Meter Info

Shows meter information.

Make support file

If you want to report a problem, you may use this.

Check update

Look if there is a later version of Ocean 2014.

Activate Meter

Use this if you want to upgrade from Connect to Professional. You will be asked to enter the license code.

Show error report

All error reports that are generated are saved in a folder on your computer. Directly when one is generated, you are always asked if you want to send it to RTI. You may not be able to do that and it will be available in this folder. Click on this button to open the folder.

About

See version information.

You can read more about Ocean 2014's program options in the topic Program options.

1.13 Close Quick Check

You can close Quick Check in two different ways:

- Close and go back to Ocean 2014's main screen
- Close and quit Ocean 2014 at the same time

Close Quick Check and go back to Ocean 2014's main screen Click on the close button:



You can also close any of the menus shown in Quick Check. This will also take you back to Ocean 2014's main screen:

Select	modality
•	Radiography
•	Fluoroscopy
•	Mammography
•	Dental

Close and quit Ocean 2014

If you want to close Ocean 2014 directly from Quick Check, click on the close button in the upper right corner of the Quick Check screen:

Preview/Print Export	Options	Close	e X Pelp
			0
110-			

Remember that Quick Check will not ask you if you want to save data unless you already have saved them once.

If you have enabled "Confirm exit" i Options, a dialogue is shown (otherwise Ocean 2104 quits directly):



You can chose to go to Windows, turn off the computer or go back to Ocean 2014.