


# DETECTOR DATA

## T20B and Piranha T20 Dose Probe



INDEPENDENT X-RAY  
QUALITY ASSURANCE

## GENERAL

The detectors **T20B Dose Probe** and **RTI T20 Dose Probe**, are all referred to as **T20** in the included documentation.

The **T20** is a type of dose detector compatible to previous detectors supplied by RTI like the R100B and the Piranha/RTI Dose Probe. It shall fulfil the relevant parts of the IEC 61674 standard for diagnostic dosimeters.

The **T20B** is optimized for the Barracuda line of electrometers with Lemo connector. The **RTI T20** is adapted for use with the Piranha/Cobia which uses a different connector.

It is optimised to contain as little material as possible, to avoid any influence on the feedback of an Image Intensifier. The long rod makes it easy to position even on the largest Image Intensifiers.

## TECHNICAL DESCRIPTION

**No corrections needed:** The detector's main part is a PIN semiconductor photo diode. The diode detects the X-rays and produces a current proportional to the intensity of the X-rays. As the diode responds differently when exposed to X-rays with different spectral content, an energy filter is placed in front of the detector to compensate for this energy dependence. This makes the response energy independent even for very low dose rates. No corrections are applied, or needs to be applied.

**EMI protected:** To make the detector more immune to EMI (ElectroMagnetic Interference), as the triaxial cable is not 100 % shielded, ferrite core(s) is used. This is placed on the end of the cable close to the connector.

**Connects to electrometer:** The currents, and charge, produced by the photo diode can then be measured with an electrometer and converted to values of the air kerma or air kerma rate (or dose and dose rate).


## SPECIFICATIONS, TYPE T20

Typical energy dependence	+1,6 to -0,4 % (radiation quality R1) +1,9 to -1,5 % (radiation quality RQR)
Air kerma (Dose)	700 pGy – 10 kGy (80 nR – 11 MR) (with Barracuda or Piranha)
Air kerma rate (Dose rate)	27 nGy/s – 500 mGy/s (3 µR/s – 57 R/s)
Inaccuracy	< 3 % between 50 – 150 kV
Typical sensitivity	8 µC/Gy (factor=125 kGy/C) (radiation quality R1)
Cable length	200 cm
Minimum cable bend radius	30 mm (∅60 mm)
Sensitive volume	This is marked by a rectangle on the red label and by a line on the side of the metallic case.
Dimensions (without cable)	35 × 320 × 6,0 mm (metal 23,5 × 5,6 mm)
Weight	Approx. 75 g
Operating temp. and humidity	-10 °C to +50 °C at <85 % relative humidity
Temperature and air pressure	Insensitive, no TP correction is necessary
Back scatter / Angular dependence	Insensitive to backscatter, will measure the same independent of what is behind the T20. Very close to ideal Cosine (flat sensor). See the graphs on page 4.

Above specifications are valid for + 18 °C to + 23 °C at <80 % relative humidity. All specifications can be altered without notice.


## OVERVIEW

The picture below shows the various parts of the T20.



## APPLICATION TIPS


The back stabilising device may be removed, as shown below, for instance when measuring in narrow places. To fit it back do the opposite.



For use of the probe with the Barracuda and the Piranha, see respective product manual. The T20 is designed for measurements in free air. For measurements behind phantoms, we suggest using the R100B or the Piranha Dose Probe. Below you can see some application photos of the T20. In the X-ray image, note the relative transparency of the rod, the black rectangular part is the metallic active detector.



T20 is taped to an Image Intensifier



X-ray image of the T20


## ENERGY CORRECTION GRAPHS

On the next page, typical correction tables and graphs for the T20 are shown, as well as graphs of the angular sensitivity.


T20		Typical Energy Correction	R1
kV	Cq	<b>Application:</b> Skin Dose, Radiographic <b>kV range:</b> 50-150 kV <b>Anode:</b> Tungsten <b>Filtration:</b> 3,0 mm Al <b>Reference point:</b> 70 kV  <b>Example:</b> Measured or set kVp = 130 kV Measured dose = 0,50 mGy Correction factor Cq = 1,002 for 130 kV Corrected dose = 1,002 x 0,50 = 0,501 mGy  <b>Inaccuracy:</b> The inaccuracy of the typical correction factor is less than ±4 %.  <b>Note:</b> Check your manual to verify when you have to do the correction manually.	
50	1,016		
55	1,010		
60	1,006		
65	1,002		
70	1,000		
75	0,998		
80	0,997		
85	0,996		
90	0,996		
95	0,996		
100	0,996		
105	0,997		
110	0,997		
115	0,998		
120	0,999		
125	1,001		
130	1,002		
135	1,004		
140	1,006		
145	1,008		
150	1,011		

T20		Typical Energy Correction	R3 (RQR)
kV	Cq	<b>Application:</b> Skin Dose, Radiographic <b>kV range:</b> 50-150 kV <b>Anode:</b> Tungsten <b>Inherent filtration:</b> 2,5 mm Al <b>Reference point:</b> 70 kV (HVL=2,58 mm Al) <b>IEC 61267 name:</b> RQR  <b>Example:</b> Measured or set kVp = 130 kV Measured dose = 0,50 mGy Correction factor Cq = 0,991 for 130 kV Corrected dose = 0,991 x 0,50 = 0,496 mGy  <b>Inaccuracy:</b> The inaccuracy of the typical correction factor is less than ±4 %.  <b>Note:</b> Check your manual to verify when you have to do the correction manually.	
50	1,019		
55	1,013		
60	1,008		
65	1,004		
70	1,000		
75	0,996		
80	0,993		
85	0,991		
90	0,988		
95	0,987		
100	0,986		
105	0,985		
110	0,985		
115	0,986		
120	0,987		
125	0,988		
130	0,990		
135	0,992		
140	0,995		
145	0,999		
150	1,002		


Typical Correction Graph



Typical Correction Graph




T20 @ 70 kV



Angular (rotational) sensitivity for the T20. The black curve is the ideal flat response.

T20 Polar Sensitivity - 70 kV



Polar plot of the same graph with the T20 located in the centre. The black curve is the ideal flat response.