



# Product Note



## Piranha – News in Ocean Next Release v4.2.1.0

### Release information

Release date	2024-09-09
Release name	Ocean Next v4.2.1.0
Firmware version	2021.9A / 4.3A

*The firmware version is available for Black and Red Piranha.*

**Note:** This firmware version does not work with old versions of Ocean 2014, or other RTI software tools, including QABrowser. See also document, *PN\_Pir\_002\_Compatibility\_Chart\_Public\_v4.2.1.0*.

### Known Limitations

None.

For further information about limitations and resolved problems, see product note *PN\_Ocean\_Next-v4.2.1.0\_Ocean*.

\*\*\* END \*\*\*