



# Product Note



## RTI Scatter Probe

News in Ocean Next Release v4.2.1.0

### Release information

Release date	2024-09-09
Release name	Ocean Next v4.2.1.0
Firmware version	v1.2
DII Version	v1.1

**Note:** This product does not work with old versions of Ocean 2014, or other RTI software tools, including QABrowser. Ocean Next v1.1.0.0 or higher is required.

This release includes various bugfixes and support for HW version 1.0 (added buzzer sound amplifier).

### Known Limitations

None.

\*\*\* END \*\*\*