



Product Note



Piranha – News in Ocean Next Release v4.2.2.0

Release information

Release date	2024-10-30
Release name	Ocean Next v4.2.2.0
Firmware version	2021.9A / 4.3A

The firmware version is available for Black and Red Piranha.

Note: This firmware version does not work with old versions of Ocean 2014, or other RTI software tools, including QABrowser. See also document, *PN_Pir_002_Compatibility_Chart_Public_v4.2.1.0*.

Known Limitations

None.

For further information about limitations and resolved problems, see product note *PN_Ocean_Next-v4.2.1.0_Ocean*.

*** END ***