



Product Note

Cobia – News in Ocean Next Release v4.4.0.0

Release information

| | |
|------------------------|---------------------|
| Release date | 202-07-10 |
| Release name | Ocean Next v4.4.0.0 |
| Cobia System version | 2024.11A |
| Cobia Firmware version | 5.20A |



New Customer Values

New firmware with support for Cobia with updated display.

Known Limitations

The table below describes limitations that have been observed and not been fixed in the release.

Description

Using RTI Chamber Adapter with multiple ion chambers, the Cobia cannot handle chambers with identical serial numbers. The result is that even if selecting the second chamber, Cobia will assume that the first one is in use. This will cause measuring error if present. However normal use of RTI chamber adapter will not include ion chambers with identical serial numbers.

*** END ***