



# Product Note

## Cobia – News in Ocean Next Release v4.6.0.0

### Release information

Release date	2026-05-05
Release name	Ocean Next v4.6.0.0
Cobia System version	2024.11A
Cobia Firmware version	5.20A



### New Customer Values

No new firmware in this release.

### Known Limitations

The table below describes limitations that have been observed and not been fixed in the release.

#### Description

Using RTI Chamber Adapter with multiple ion chambers, the Cobia cannot handle chambers with identical serial numbers. The result is that even if selecting the second chamber, Cobia will assume that the first one is in use. This will cause measuring error if present. However normal use of RTI chamber adapter will not include ion chambers with identical serial numbers.

\*\*\* END \*\*\*