

Product Note

Mako – News in Ocean Next Release v4.6.0.0



Release information

Release date	2026-05-05
Release name	Ocean Next v4.6.0.0
Firmware version	1.6.4
Mako DLL	1.5.291

Note: This firmware version does not work with old versions of Ocean 2014, or other RTI software tools, including QABrowser. See also document, *PN_Pir_002_Compatibility_Chart_Public_v4.6.0.0*.

Mako, the most accurate, efficient and versatile X-ray meter ever. The next-generation multimodality X-ray meter has market-leading accuracy ($\pm 1.5\%$ kVp uncertainty) thanks to innovative new detector technology. Centered around efficiency, Mako provides the ultimate no-fuss experience, simple setup in all X-ray applications, and fully wireless capability.

The modular, future-proof system design provides ultimate configuration flexibility, with the widest-ever application range for X-ray Quality Assurance testing.

This release fixes several bugs.

Known Limitations

None.

For further information about limitations and resolved problems, see product note *PN_Ocean_Next-v4.6.0.0_Ocean*.

*** END ***