



# Product Note

## Cobia – News in Release R16-2.7.1

### Release information

Release date	2020-02-19
Release name	R16-2.7.1
Cobia System version	2019.12A
Cobia Firmware version	5.13C



### New Customer Values

No new Cobia firmware with this release.

### Known Limitations

The table below describes limitations that have been observed and not been fixed in the release.

#### Description

Using RTI Chamber Adapter with multiple ion chambers, the Cobia cannot handle chambers with identical serial numbers. The result is that even if selecting the second chamber, Cobia will assume that the first one is in use. This will cause measuring error if present. However normal use of RTI chamber adapter will not include ion chambers with identical serial numbers.

\*\*\* END \*\*\*