



Product Note



Piranha – News in Release R16-2.7.1

Release information

Release date	2020-02-19
Release name	R16-2.7.1
Installer version	2020.2A
Firmware version	4.2D

The firmware version is available for both Black and Red Piranha. However, some functionality available in Black Piranha is not available in Red Piranha, such as new radiation qualities.

This software release for Piranha comes together with releases of Ocean 2014 for Cobia and Black Piranha.

Note: This firmware version does not work with old versions of Ocean 2014, or other RTI software tools, including QABrowser. See also document, *PN_Pir_002_Compatibility_Chart_Public_v2020.2A*.

New Customer Values

No new changes in this release.

Known Limitations

None.

For further information about limitations and resolved problems, see product note *PN_R16-2.7.1_Ocean*.

*** END ***