




INDEPENDENT X-RAY
QUALITY ASSURANCE

Use with

Piranha Cobia



Ocean Next™ general specifications

The new power in X-ray QA software

Since its launch, Ocean™ has been *the* market-leading software for X-ray QA, addressing customers' needs all the way from Quick, for daily checks, to advanced QA with workflow, automation, and traceability in Professional in easy steps.

If you want a swift, easy-to-use application for routine controls, or a customized application with workflow, automatic tests, and traceability, the new Ocean Next™ software suite provides you with both options.

Its recognition and efficiency make it central to our X-ray solutions, and it is compliant with all our QA meters.

You can plan the measurements in advance, create checklists, add information as a pop-up window for a specific exposure, and include instructions to simplify the work for you and your co-workers using a streamlined user interface.

Choose from three different license levels: *Quick*, *Advantage*, and *Professional*.

Key features

- Customized application with workflow, automatic tests, and traceability
- Updated user interface to recognition and efficiency
- Enhanced functionality in Quick Check
- Upgraded with simplified templates and measurements structure
- Runs on both computers and tablets - choose the display you want!

RTI Group Headquarters
Flöjelbergsgatan 8C
SE-431 37 Mölndal, SWEDEN

Phone: +46 (0) 31 746 36 00
info@rtigroup.com
www.rtigroup.com

Which Ocean Next™ suits you best?

Built-in display function	Quick/Advantage/Professional
Save function	Quick/Advantage/Professional
Basic report printing	Quick/Advantage/Professional
Basic Excel connection	Quick/Advantage/Professional
Advanced report printing	Advantage/Professional
Advanced Excel connection	Advantage/Professional
User-defined tests & checklists	Advantage/Professional
Analysis with pass/fail limits	Advantage/Professional
Reporting	Advantage/Professional
CT Dose Profiler support	Advantage/Professional
Single-page templates	Advantage/Professional
Multi-page templates	Professional
Site database	Professional

<https://rtigroup.com/product-lines/software/>