



# Product Note

## Ocean Next – News in Release v1.0.0.0

### Release information

Release date	2021-01-11
Release name	v1.0.0.0
EXE file version	1.0.2.26165
DLL version	3.0.1.9132

This is a Major Feature release of Ocean Next.

Major changes includes updated graphical user interface, multiple documents support, new license levels, improved scaling for different DPIs and resolutions, bug fixes, improved performance and many more things. Refer to <https://rtigroup.com/product-lines/software/> for a full list of features in the new Ocean Next software.

Note that this version of Ocean installs per user. Ocean Next will only be available to the current windows user that was logged in when installing Ocean.

### Database format

This version of Ocean has a new database structure compared to Ocean 2014.

Note that the location of the Ocean database has moved to:  
"C:\Users\[WindowsUser]\AppData\Local\RTI Group\Ocean Next\ProgramData\AllData.accdb"

### Cobia and Piranha firmware

#### Piranha

No new Piranha firmware version, in this release. The lowest firmware that shall be used with this release of Ocean is 4.2D.

#### Cobia

No new Cobia firmware version, in this release. The lowest firmware that that is recommended to be used with this version of Ocean Next is 2019.12A.



# Product Note

## Known Limitations

The table below shows issues that have been verified, and valid in this release of Ocean Next.

Description
Under rare circumstances the Piranha may lose the Ocean license. Contact support (support@rtigroup.com) if this occurs to resolve issue.

\*\*\* END \*\*\*