



# Product Note



## Piranha – News in Ocean Next Release v2.0.0.1

### Release information

Release date	2022-05-02
Release name	Ocean Next v2.0.0.1
Firmware version	2021.9A / 4.3A

*The firmware version is available for Black and Red Piranha.*

Major changes include:

- Added support for new Bluetooth module.
- Updated Dose Correction for RTI Dose Probe for Mammography Applications
- Renamed some Beam Qualities
- New CT Calibrations
- New Mammo Calibrations
- Added protection for memory corruption

**Note:** This firmware version does not work with old versions of Ocean 2014, or other RTI software tools, including QABrowser. See also document, *PN\_Pir\_002\_Compatibility\_Chart\_Public\_v2.0.0.1.*

## Known Limitations

None.

For further information about limitations and resolved problems, see product note *PN\_Ocean\_Next-v2.0.0.0\_Ocean.*

\*\*\* END \*\*\*