



Product Note



Piranha – News in Ocean Next Release v2.0.1.0

Release information

Release date	2022-05-13
Release name	Ocean Next v2.0.1.0
Firmware version	2021.9A / 4.3A

The firmware version is available for Black and Red Piranha.

Major changes include:

- Added support for new Bluetooth module.
- Updated Dose Correction for RTI Dose Probe for Mammography Applications
- Renamed some Beam Qualities
- New CT Calibrations
- New Mammo Calibrations
- Added protection for memory corruption

Note: This firmware version does not work with old versions of Ocean 2014, or other RTI software tools, including QABrowser. See also document, *PN_Pir_002_Compatibility_Chart_Public_v2.0.0.1.*

Known Limitations

None.

For further information about limitations and resolved problems, see product note *PN_Ocean_Next-v2.0.0.0_Ocean.*

*** END ***