



# Piranha Upgrade

SER-011\_bil.3\_A



**Piranha S/N CB2-** .....

## Company Information

Your Company Name:
Contact person:
E-mail:
Phone:
Order No:

## Invoice Address

## Shipping Address *If different from Invoice Address*

Company name:		Company name:	
Contact person:		Department:	
Street:		Street:	
Zip/post code:	City:	Zip/post code:	City:
Country:		Country:	

End user information on certificates	Standard Freight <i>(RTI handles freight and adds charges to invoice.)</i>
Company Name:	
Street:	<b>Own Forwarder:</b>
Zip/post code:	Freighter name:
City:	Freight Account No:
Country:	

The RTI Calibration Laboratory in Mölndal, Sweden, is an Accredited Dosimetry Calibration Laboratory.

The calibrations are traceable to international standards and accredited by SWEDAC to ISO/IEC 17025:2017. SWEDAC is a member of the EA, IAF, and ILAC.

# Piranha Upgrade options

## Piranha Hardware Upgrade:

### 9829302-02 Red to Black (R2B)

#### Upgrade your Red Piranha to a Black Piranha!

Improved Bluetooth connectivity with auto-search and pairing between Ocean and Piranha.

Connect to your tablet or PC via built-in Bluetooth or USB. The Black Piranha is complemented by the Ocean Next database software. Black Piranha automatically recognizes probes being used and selects appropriate settings for those measurements. Ocean Next also features Quick Check mode, a quick and easy way to start your measurements. Simple and accurate measurements – data and waveforms are instantly displayed. Measurement templates are provided, just power-up and go.

Black Piranha offers an improved, leadfree, solid- state detector design provides for greater accuracy and ease in positioning.

#### Included:

New Piranha Built-in battery.

Calibration cost for the main unit.

Extended Warranty Program 24 months on the Piranha main unit

*(Possible to Prolonged extended warranty not later than 24 months from last calibration.)*

Ocean Next Advantage license level is included, unless the Piranha already has Professional license.

The Piranha retains its serial number *(but is essentially rebuilt to a new instrument!)*.

## Piranha Modality Upgrade:

### 9620557-00 Piranha Modality Upgrade from Piranha D 2xx (Dental) to Piranha R/F (557)

Add Radiography, Fluoroscopy, and CT to your Piranha.

Upgrade any Dental Piranha model (Piranha D 200, 201, 204, 205, 250, 251, 254, 255, Piranha Dental) to a Piranha R/F (used to be called Piranha 557).

### 9620657-00 Piranha Modality Upgrade from Piranha 2xx (Dental) to Piranha MULTI (657)

Add Radiography, Fluoroscopy, Mammography, and CT to your Piranha.

Upgrade any Dental Piranha model (Piranha D 200, 201, 204, 205, 250, 251, 254, 255, Piranha Dental) to a Piranha Multi (used to be called Piranha 657).

### 9630657-00 Piranha Modality Upgrade from Piranha 3xx (Mammo) to Piranha MULTI (657)

Add Radiography, Fluoroscopy, CT, and Dental to your Piranha.

Upgrade any Mammo Piranha model (Piranha M 300, 301, 304, 305, 350, 351, 354, 355, Piranha Mammo) to a Piranha Multi (used to be called Piranha 657).

### 9640557-00 Piranha Modality Upgrade from Piranha 4xx (CT) to Piranha R/F (557)

Add Radiography, Fluoroscopy, and Dental to your Piranha.

Upgrade any CT Piranha model (Piranha CT 400, 401, 450, 451, Piranha CT) to a Piranha R/F (used to be called Piranha 557).

### 9640657-00 Piranha Modality Upgrade from Piranha 4xx (CT) to Piranha Multi (657)

Add Radiography, Fluoroscopy, Mammography, and Dental to your Piranha.

Upgrade any CT Piranha model (Piranha CT 400, 401, 450, 451, Piranha CT) to a Piranha Multi (used to be called Piranha 657).

### 9650657-00 Piranha Modality Upgrade from Piranha 5xx (R/F) to Piranha MULTI (657)

Add Mammography to your Piranha.

Upgrade any R/F Piranha model (Piranha R&F 500, 501, 504, 505, 506, 507, 550, 551, 554, 555, 556, 557, Piranha R/F) to a Piranha Multi (used to be called Piranha 657).