




INDEPENDENT X-RAY
QUALITY ASSURANCE

Use with

- Mako RTI Scatter Probe Piranha Cobia



Ocean Next™ version 4.0 General specifications

The world-leading software for X-ray QA & testing

Since its launch in 2010, Ocean™ has been the world-leading software for X-ray QA and testing, addressing customers' needs from Quick, for daily checks, to the advanced QA with workflow, automation, and traceability in Professional – in easy steps.

The latest version, Ocean Next software* v.4.0, introduces upgraded and new automated features, integration with the new Mako, myRTI, and more. Whether you require a swift, user-friendly application for routine controls or a customized application with workflow, automatic tests, and traceability, Ocean offers both options. Ocean Next 4.0 also brings major enhancements, including a smoother workflow with faster response time, new product support features, and digital services.

Ocean plays a central role in our X-ray solutions, ensuring compliance with our multimodality Mako, Piranha and Cobia meters, as well as the RTI Scatter Probe for combined leakage and scatter detection. With a streamlined user interface, you can plan measurements in advance, create checklists, add pop-up information for specific exposures, and include instructions to simplify work for you and your co-workers.

Choose from three different license levels Quick, Advantage, or Professional. Visit [Ocean Next software](#) on [rtigroup.com](#).

Don't forget to register for a myRTI account!

*Ocean Next requires MS Windows 10 or later.

Key features

- Customized application with workflow, automatic tests, and traceability.
- Updated, streamlined user interface for recognition and efficiency.
- Message panel for RTI news and more visible, automated calibration reminders.
- Calibration notices available in reports.
- Simplified templates and measurements structure.
- Runs on both computers and tablets – choose the display you want.
- RTI Academy training courses available.

Which Ocean Next™ suits you best?

Built-in display function	Quick/Advantage/Professional
Save function	Quick/Advantage/Professional
Printing	Quick/Advantage/Professional
Basic Excel® connection	Quick/Advantage/Professional
Advanced Excel® connection	Advantage/Professional
User-defined tests & checklists	Advantage/Professional
Analysis with pass/fail limits	Advantage/Professional
Reporting	Advantage/Professional
CT Dose Profiler support	Advantage/Professional
Single-page templates	Advantage/Professional
Multi-page templates	Professional
Site database (incl. location, room, equipment, etc.)	Professional
Scheduling	Professional



RTI Group Headquarters
Flöjelbergsgatan 8C
SE-431 37 Mölndal, SWEDEN

Phone: +46 (0) 31 746 36 00
sales@rtigroup.com
www.rtigroup.com

RTI Group North America
33 Jacksonville Road, Bldg. 1
Towaco, NJ 07082, USA

Phone: +1 800-222-7537
sales.us@rtigroup.com
www.rtigroup.com