



INDEPENDENT X-RAY
QUALITY ASSURANCE

Use with



Mako Suction Cup Holder | Specifications

Simple Positioning of Mako Probes on Generic Surfaces

The RTI Mako Suction Cup Holder allows for easy positioning of the Mako Probes (R/F, Mammo or Dental) on generic surfaces, such as detector plates or covers, in vertical or horizontal positions.

Application

The Mako Probes can be easily inserted to the Mako Suction Cup Holder and connected to the Mako Base Unit via USB cable. The Suction Cups can then be used for firm adhesion to a flat surface, to allow flexible and simple positioning of the Mako Probe in the X-ray path.

Caution: ensure firm adhesion to the required surface to avoid dropping and/or damaging the Mako Probe.



Specifications

| | |
|--------------------|-------------------------|
| Article number: | 9765038-00 |
| Article name: | Mako Suction Cup Holder |
| Dimensions: | 107 x 40 x 21,5 mm |
| Weight: | 30 g |
| Material: | POM |
| Compatible meters: | Mako Probes |

